

Back Wrap = Sleeve + 2 Heat Exchangers + Y Connector

Your Equine Back Wrap has been packaged with each of the two removable Heat Exchangers already inserted into the Sleeve. Before you remove the Heat Exchangers for cleaning, take care to note their assembly in the Sleeve. The Y-Connector is bagged separately.

Fig. 1 Back Wrap Heat Exchanger Removal

1. Disconnect the Y Connector from the Heat Exchangers.
2. Place the Sleeve on a flat surface with the logo side down. Fold back the flaps along the center of each side of the sleeve to open the zippers. Gently pull out each Heat Exchanger from each side of the Wrap.







Fig. 2 Back Wrap Heat Exchanger Insertion

1. Place the Sleeve on a flat surface with the logo side down. Fold back the flaps along the center of each side of the sleeve to open the zippers. Insert each of the two Heat Exchangers into the corresponding side of the Wrap with the black side facing down.
2. Slide each Heat Exchanger Connector through the corresponding opening in the Sleeve. The part of the Sleeve which you have folded over to open the zipper should be pulled over the top of the Heat Exchanger to close the zipper. The green side of the Heat Exchanger should touch the green side of the Sleeve.



Important: Make sure the Heat Exchanger is flat inside the Sleeve with no folds or creases. After assembly, fold at the natural horizontal fold, and make sure the Wrap is flat and smooth.




IMPORTANT SAFETY INFORMATION

- 1 Always consult with your veterinarian before using the Game Ready Equine System for an injury.
- 2 During treatment, the horse must be accompanied by an attendant.
- 3 Do not place the Wrap directly against open wounds, sores, rashes, infections, or stitches. First place a sterile dressing on the wound and then apply the Wrap.
- 4 Do not store in a closed plastic bag after use, as this may promote mold growth.
- 5 Do not puncture the Wrap; doing so will void the Game Ready Equine warranty.
- 6 Initial applications of the Back Wrap may cause some horses mild concern at the new sensation. Carefully monitor your horse's reactions throughout each step of the process and proceed only when your horse appears calm and relaxed.
- 7 It is **recommended** that the highest level of pressure is Medium Pressure (~~~~) for this Wrap. If the pressure appears uncomfortable, pause the Unit, and restart with no pressure.

CoolSystems, Inc. 1201 Marina Village Pkwy, Suite 200 Alameda, CA 94501 1-888-GAMEREADY
www.gamereadyequine.com

Back Wrap Application

Note: For maximum cooling, place the Wrap against the horse's bare back, unless there is a wound. In the case of a wound, cover with a bandage before applying the Wrap.



Step 1: Place the Wrap on the horse's back, aligning the hose connection on top of the horse's withers, similar to the positioning of a pommel or saddle horn. Draw the front strap from the far side of the horse, under his belly, clipping to the nearside front buckle.

Step 2: Adjust the tautness by shortening the nylon strap through the sliding buckles on both near and far sides. Strap should be snug to the horse's belly, but not as tight as a girth. Follow the same protocol for the rear strap, carefully monitoring your horse's reaction. Do not proceed further until your horse appears calm and relaxed.

Step 3: Clip the upper Breastplate Straps into either side of the Back Wrap at the horse's shoulders. Pass the bottom Breastplate Strap through the horse's front legs and clip to the Front Girth Strap. Adjust the tautness by shortening the nylon straps at all points.


Step 4: Use the lacing along the spine of the Back Wrap to customize its shape for optimal fit and contact. For best contact, adjust each side of the Back Wrap to cover the horse's spine. The right and left sides of the Wrap may overlap slightly, but should not be folded over onto themselves.

Step 5: Clip the split end of the Y Connector to each of the two Heat Exchanger Connectors. Clip the single end of the Y Connector to the Single Hose Connector. Clip the Single Hose Connector to the Control Unit.

Begin initial treatment with no pressure in order to acclimate the horse to the sensation of the Back Wrap. If after several minutes there are no signs of aversion from the horse, begin administering cold therapy at low pressure. Advance to medium pressure only when your horse appears comfortable and relaxed. **[Medium pressure should be the highest setting used for the Back Wrap, as the product performs optimally at this setting.]**

Caution: First application of this Wrap may cause horses to startle preliminarily at the new sensation. Proceed with administering therapy with caution. Do not open the zipper on the Wrap while the Wrap is in use, or while being applied to the horse.

Protocol: 30-40 minute application on low, medium, or no pressure. Always consult your veterinarian before using the Game Ready Equine System on any injury.



Note:

Run the Connector Hose through the Strap to fix the Connector in place.

Do not allow the Hose to bend the Wrap at the connection. This can restrict water flow and affect cooling.

After treatment, disconnect the Hose from the Wrap first, before removal of the Wrap from the horse.

Care: Gently remove the Heat Exchangers from the Sleeve. Hand wash the Sleeve in cold water and mild detergent, or antibacterial soap. Hang to dry. Hand wash the Heat Exchanger, do not machine wash; doing so will void the warranty. Periodically run hydrogen peroxide through your Game Ready Equine System to keep the interior water chamber of your Wrap clean. See Control Unit User's Manual for more detail. Cold sterilization only.

Storage: Store the Back Wrap laying flat with green sides out or hang. Do not bend lengthwise.

Warranty Information:

Sleeve: In case of manufacturer defect: Sleeve may be returned within 7 days of purchase.

Heat Exchanger & Y-Connector: : 1 year from date of purchase. See warranty card included .